Himalayan SALT PAD

For Hot & Cold Therapy

La Naturess heating pad is designed with high-quality natural cotton and filled with 100% Himalayan Salt crystals. Salt Heat therapy is in use as a traditional natural remedy for pain relief since ancient times. Himalayan salt has heat absorption properties which help to retain the heat in salt and keep the heating pad warm for a longer time. The Salt pad can also be used for cold therapy due to its low-temperature point and holds any cold or hot temperatures. Heat in the microwave for heat therapy or chill in the freezer for cold therapy.

For Microwave Use:

- Place the heat pad in microwave
- Under Medium setting, heat for 2 min and add increments of 20 seconds for extra heat.

For Cold Therapy:

- To use as a cold pack put the pad inside a sealable plastic bag in the freezer for 1-2 hours.

Benefits:

- Relieve muscle spasms, pain, inflammation, and edema
- Decrease joint stiffness, menstrual cramps, and arthritis pain
- Increase extensibility of collagen tissues
- Increase blood flow and help with an earache, toothache and muscle pain

Precautions:

- Follow the Heating Instructions and DO NOT overheat product
- Clean microwave from any previous spills and oil which can cause the product to burn
- Do NOT heat the product for more than 3 min continuously
- Check the heat pad temperature before applying it to the body
- DO NOT wash the head pad with the white salt bag inside it. Remove the white bag and ONLY wash the removable cotton bag in Green.
- Dry the heating pad's cover carefully before using it again.
- Keep out of reach of children
- Keep Product away from water

Contact Information: Coordonnées

38- 1200 Aerowood Drive
Mississauga,
ON. Canada. L4W 2S7

905-625-6530
905-510-1165

Email: info@lanaturess.com / Web: www.lanaturess.com

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